

TAURUSDISPO® Cell Culture Plate provides uniform and compatible surfaces for animal cell attachment and growth. TAURUSDISPO® Cell Culture Plates suit a multitude of experimental needs while providing cost efficiency for a wide selection of applications in your cell culture routine. TAURUSDISPO® Cell Culture Plate is individually wrapped with clear lot number for batch traceability. The TC Treated Cell Culture Plate has crystal-clear surfaces for distortion-free microscopy and multiple features that help reduce and prevent cross-contamination




**UNIFORM
HOLE POSITION**



**100%
VIRGIN-POLYSTYRENE**



**STERILIZED
BY E-BEAM**



**INDIVIDUAL
PACKAGING**

- High clarity, 100% virgin polystyrene.
- "TC" refers to vacuum plasma tissue culture treatment.
- Available for 6, 12, 24, 48 and 96 well plates.
- Sterilized by irradiation, SAL=10⁻⁶.
- Individually Wrapped and Lot No. Traceable



Cat.No.	Well	Description	Bottom Type	Color	Material	Packaging
TI-CMP-006	06-Well	TC Treated, Individually Wrapped, Sterile	Flat Bottom	Clear	PS	50 packs/case
TI-CMP-012	12-Well	TC Treated, Individually Wrapped, Sterile	Flat Bottom	Clear	PS	50 packs/case
TI-CMP-024	24-Well	TC Treated, Individually Wrapped, Sterile	Flat Bottom	Clear	PS	50 packs/case
TI-CMP-048	48-Well	TC Treated, Individually Wrapped, Sterile	Flat Bottom	Clear	PS	50 packs/case
TI-CMP-096	96-Well	TC Treated, Individually Wrapped, Sterile	Flat Bottom	Clear	PS	50 packs/case
TI-CMP-096U	96-Well	TC Treated, Individually Wrapped, Sterile	U - Bottom	Clear	PS	50 packs/case
TI-CMP-106	06-Well	Non Treated, Individually Wrapped, Sterile	Flat Bottom	Clear	PS	50 packs/case
TI-CMP-112	12-Well	Non Treated, Individually Wrapped, Sterile	Flat Bottom	Clear	PS	50 packs/case
TI-CMP-124	24-Well	Non Treated, Individually Wrapped, Sterile	Flat Bottom	Clear	PS	50 packs/case
TI-CMP-148	48-Well	Non Treated, Individually Wrapped, Sterile	Flat Bottom	Clear	PS	50 packs/case
TI-CMP-196	96-Well	Non Treated, Individually Wrapped, Sterile	Flat Bottom	Clear	PS	50 packs/case
TI-CMP-196U	96-Well	Non Treated, Individually Wrapped, Sterile	U - Bottom	Clear	PS	50 packs/case

Kindly Reminder:

Please feel free to contact us.
For further information or sample ordering please visit our website
www.taurusinc.co.in or contact us: sales@taurusinc.co.in